

Protocol for EQAsia EQA4 2022

ID and antimicrobial susceptibility testing of *Klebsiella pneumoniae*, *Acinetobacter* spp. and *Staphylococcus aureus* test strains

1	INTRODUCTION	1
2	OBJECTIVES	2
3	OUTLINE OF THE EQASIA EQA	2
3.1	Shipping, receipt and storage of strains	2
3.2	Identification of <i>Klebsiella pneumoniae</i> , <i>Acinetobacter</i> spp. and <i>Staphylococcus aureus</i> test strains	3
3.3	Antimicrobial susceptibility testing of <i>Klebsiella pneumoniae</i> , <i>Acinetobacter</i> spp. and <i>Staphylococcus aureus</i> test strains, and of the reference strains	3
4	REPORTING OF RESULTS AND EVALUATION	9
5	HOW TO SUBMIT RESULTS VIA THE INFORMATICS MODULE	9

Version 2 of the protocol includes changes in Table 1 reference values (marked in bold).

1 INTRODUCTION

The EQAsia project aims to strengthen the provision of External Quality Assessment (EQA) services across the One Health sector in South and Southeast Asia. Therefore, a comprehensive and high-quality EQA program for antimicrobial resistance (AMR) is offered to all the National Reference Laboratories/Centres of Excellence in the region during 2021. The EQA is organized by the consortium of EQAsia and supported by the Fleming Fund.

The EQAsia EQA4 2022 includes the antimicrobial susceptibility testing of five *Klebsiella pneumoniae*, five *Acinetobacter* spp. and five *Staphylococcus aureus* strains **identified** among a total of **seven** test strains for each microorganism, which include two non-target species strains.

Additionally, antimicrobial susceptibility testing of the relevant reference strains for quality control (QC) in relation to antimicrobial susceptibility testing is included. The QC reference strains supplied (or that have been supplied in previous EQAS) are: *Escherichia coli* ATCC 25922/CCM 3954, *E. coli* NCTC 13846/CCM 8874 (for colistin), *Pseudomonas aeruginosa* ATCC 27853/CCM 3955, *S. aureus* ATCC 25923/CCM 3953 (for disk diffusion) and *S. aureus* ATCC 29213/CCM 4223 (for MIC). These reference strains are original CERTIFIED cultures provided free of charge, and should

be used for future internal quality control for antimicrobial susceptibility testing in your laboratory. Therefore, please take proper care of these strains. Handle and maintain them as suggested in the manual 'Subculture and maintenance of quality control strains' available on the [EQAsia website](#).

2 OBJECTIVES

The main objective of this EQA is to support laboratories to assess and if necessary, improve the identification and antimicrobial susceptibility testing of pathogens, specifically *Klebsiella pneumoniae*, *Acinetobacter* spp. and *Staphylococcus aureus*. Therefore, the laboratory work for this EQA should be performed using the methods routinely used in your own laboratory.

3 OUTLINE OF THE EQASIA EQA

3.1 Shipping, receipt and storage of strains

In April 2022, it is expected that around 25 laboratories located in South and Southeast Asia will receive a parcel containing one or more of the following:

- Seven test strains of which five are *Klebsiella pneumoniae*, in addition to two non-target species strains. The *Escherichia coli* ATCC 25922/CCM 3954 and *E. coli* NCTC 13846/CCM 8874 (for colistin) will be provided as reference strains (if not already received in previous EQAs).
- Seven test strains of which five are *Acinetobacter* spp., in addition to two non-target species strains. The *Pseudomonas aeruginosa* ATCC 27853/CCM 3955 will be provided as reference strain (if not already received in previous EQAs).
- Seven test strains of which five are *Staphylococcus aureus*, in addition to two non-target species strains. The *S. aureus* ATCC 25923/CCM 3953 (for disk diffusion) and *S. aureus* ATCC 29213/CCM 4223 (for MIC) will be provided as reference strains (if not already received in previous EQAs).

Please confirm receipt of the parcel through the confirmation form enclosed in the shipment

All strains are shipped lyophilized. The lyophilized strains must be stored in a dark, cool place. The strains must be sub-cultured and prepared for storage in your strain collection (e.g. in a -80°C freezer). This set of cultures should serve as reference if discrepancies are detected during the testing (e.g. they can be used to detect errors such as mislabelling or contamination), and they can function as reference material available for reference at a later stage, when needed.

For reconstitution of the test strains, please see the document 'Instructions for opening and reviving lyophilised cultures of test strains' on the [EQAsia website](#).

For reconstitution of the QC reference strains, please see the document 'Subculture and maintenance of quality control strains' on the [EQAsia website](#).

All provided strains belong to UN3373, Biological substance category B. These strains can potentially be harmful to humans and pose a risk due to their possible pan-resistant profile, therefore becoming

a challenge in the treatment of a potential human infection. It is the recipient laboratory's responsibility to comply with national legislation, rules and regulations regarding the correct use and handling of the provided test strains, and to possess the proper equipment and protocols to handle these strains. Nevertheless, it is recommended to handle the strains in a BSL2 containment facility using equipment and operational practices for work involving infectious or potentially infectious materials. The containment and operational requirements may vary with the species, subspecies, and/or strains, thus, please take the necessary precautions.

Please consult the [Pathogen Safety Data Sheets](#) (PSDSs) produced by the Public Health Agency of Canada. The PSDSs of each pathogen can be found in the bottom of the page. These PSDSs are technical documents that describe the hazardous properties of human pathogens, and provide recommendations for the work involving these agents in a laboratory setting.

3.2 Identification of *Klebsiella pneumoniae*, *Acinetobacter* spp. and *Staphylococcus aureus* test strains

For each test species, two out of the seven test strains related to each bacterial species does not belong to the target species of the EQA trial. To identify the five cultures of the correct target species among the seven test strains, you should use the method routinely used in your own laboratory for **identification** of the organism.

3.3 Antimicrobial susceptibility testing of *Klebsiella pneumoniae*, *Acinetobacter* spp. and *Staphylococcus aureus* test strains, and of the reference strains

The strains identified as *Klebsiella pneumoniae*, *Acinetobacter* spp. and *Staphylococcus aureus*, as well as the appropriate reference strains, should be tested for susceptibility towards as many as possible of the antimicrobials mentioned in the test form and in **Tables 1-3**. Note that some of the antimicrobials (highlighted) could be omitted by the Human Health laboratories. Please use the methods routinely used in your own laboratory.

The reference values used in this EQA for interpreting MIC and disk diffusion results are in accordance with current zone diameter and MIC breakpoint values developed by CLSI (M100, 31st Ed.). When not available, EUCAST clinical breakpoints (Tables v. 11.0, 2021) or epidemiological cut off values (<https://mic.eucast.org/>) are used instead. The breakpoint values for *Klebsiella pneumoniae*, *Acinetobacter* spp. and *Staphylococcus aureus* can be found in **Tables 1-3**, respectively. **Make sure to use the correct table for the interpretation.**

Interpretation of MIC or disk diffusion results will lead to categorization of the result into one of the categories: **resistant** (R), **intermediate** (I) or **susceptible** (S). In the evaluation report you receive upon the submission deadline, the obtained interpretation in comparison with the expected interpretation will be evaluated/scored as follows:

SCORES		Obtained Interpretation		
		Susceptible	Intermediate	Resistant
Expected Interpretation	Susceptible	4	3	1
	Intermediate	3	4	3
	Resistant	0	3	4

0	Incorrect: very major
1	Incorrect: major
3	Incorrect: minor
4	Correct

Table 1. Interpretive criteria for *Klebsiella pneumoniae* antimicrobial susceptibility testing

The highlighted antimicrobials could be omitted by the Human Health laboratories.

Antimicrobials	Reference values			Reference values		
	MIC ($\mu\text{g/mL}$)			Disk diffusion (mm)		
	S	I	R	S	I	R
Amikacin, AMK	≤ 16	32	≥ 64	≥ 17	15-16	≤ 14
Ampicillin, AMP	≤ 8	16	≥ 32	≥ 17	14-16	≤ 13
Azithromycin, AZI	≤ 16	-	≥ 32	≥ 13	-	≤ 12
Cefepime, FEP	≤ 2	4-8	≥ 16	≥ 25	19-24	≤ 18
Cefotaxime, FOT	≤ 1	2	≥ 4	≥ 26	23-25	≤ 22
Cefotaxime/clavulanic acid, F/C	NA	NA	NA	NA	NA	NA
Cefoxitin, FOX	≤ 8	16	≥ 32	≥ 18	15-17	≤ 14
Ceftazidime, TAZ	≤ 4	8	≥ 16	≥ 21	18-20	≤ 17
Ceftazidime/clavulanic acid, T/C	NA	NA	NA	NA	NA	NA
Chloramphenicol, CHL	≤ 8	16	≥ 32	≥ 18	13-17	≤ 12
Ciprofloxacin, CIP	≤ 0.25	0.5	≥ 1	≥ 26	22-25	≤ 21
Colistin, COL	-	≤ 2	≥ 4	NA	NA	NA
Doripenem, DOR	≤ 1	2	≥ 4	≥ 23	20-22	≤ 19
Ertapenem, ETP	≤ 0.5	1	≥ 2	≥ 22	19-21	≤ 18
Gentamicin, GEN	≤ 4	8	≥ 16	≥ 15	13-14	≤ 12
Imipenem, IMI	≤ 1	2	≥ 4	≥ 23	20-22	≤ 19
Levofloxacin, LEVO	≤ 0.5	1	≥ 2	≥ 21	17-20	≤ 16
Meropenem, MERO	≤ 1	2	≥ 4	≥ 23	20-22	≤ 19
Nalidixic acid, NAL	≤ 16	-	≥ 32	≥ 19	14-18	≤ 13

Piperacillin/tazobactam, PT4	$\leq 8/4$	16/4	$\geq 32/4$	≥ 25	21-24	≤ 20
Sulfamethoxazole, SMX	≤ 256	-	≥ 512	≥ 17	13-16	≤ 12
Tetracycline, TET	≤ 4	8	≥ 16	≥ 15	12-14	≤ 11
Tigecycline, TGC*	≤ 2	-	≥ 4	NA	NA	NA
Tobramycin, TOB	≤ 4	8	≥ 16	≥ 15	13-14	≤ 12
Trimethoprim, TMP	≤ 8	-	≥ 16	≥ 16	11-15	≤ 10
Trimethoprim/sulfamethoxazole, SXT	$\leq 2/38$	-	$\geq 4/76$	≥ 16	11-15	≤ 10

Reference values are based on Enterobacterales breakpoints from CLSI M100, 32nd Ed.

*Reference values are based on *K. pneumoniae* epidemiological cut off values from <https://mic.eucast.org/> on January 2022.

Beta-lactam and carbapenem resistance

The following tests for detection of ESBL-, AmpC-, and carbapenemase-producing phenotypes for *K. pneumoniae* are recommended:

- Reduced susceptibility to cefotaxime (FOT) and/or ceftazidime (TAZ): it indicates that the bacterial strain is an ESBL-, AmpC, or carbapenemase-producing phenotype. These strains should be tested for ESBL-, AmpC, or carbapenemase-production by confirmatory tests.
- Confirmatory test for ESBL production: it requires the use of both cefotaxime (FOT) and ceftazidime (TAZ) alone, as well as in combination with a β -lactamase inhibitor (clavulanic acid). Synergy can be determined by broth microdilution methods, Gradient Test or Disk Diffusion:
 - It is defined as a ≥ 3 twofold concentration decrease in an MIC for either antimicrobial agent tested in combination with clavulanic acid vs. its MIC when tested alone (Gradient Test 3 dilution steps difference; MIC FOT : FOT/Cl or TAZ : TAZ/Cl ratio ≥ 8).
 - A positive synergy testing for Disk Diffusion is defined as ≥ 5 mm increase of diameter of FOT or TAZ in combination with clavulanic acid (FOT/Cl or TAZ/Cl) compared to testing them alone. The presence of synergy indicates ESBL production.
- Detection of AmpC-type beta-lactamases: it can be performed by testing the bacterial culture for susceptibility to ceftaxitin (FOX). Resistance to FOX indicates the presence of an AmpC-type beta-lactamase.
- Confirmatory test for carbapenemase production: it requires the testing of meropenem (MERO). Resistance to MERO indicates that the bacterial strain is a carbapenemase-producer.

It should be noted that some resistance mechanisms do not always confer clinical resistance. Therefore, the classification of the phenotypic results (**Figure 1** below) should be based on the “EUCAST guidelines for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance”, Version 2.0, July 2017, and the most recent EFSA recommendations – The European Union summary report on antimicrobial resistance in zoonotic and indicator bacteria from humans, animals and food in 2017/2018, EFSA Journal 2020;18 (3) <https://doi.org/10.2903/j.efsa.2020.6007>

1. ESBL-Phenotype			4. Carbapenemase-Phenotype		
	MIC (mg/L)	Zone Diameter (mm)		MIC (mg/L)	Zone Diameter (mm)
FOT or TAZ	> 1	< 21 (FOT); < 22 (TAZ)	MERO	> 0.12	< 25
MERO	≤ 0.12	≥ 25	5. Other Phenotypes		
FOX	≤ 8	≥ 19			
FOT/CLV and/or TAZ/CLV	SYNERGY	SYNERGY			
2. AmpC-Phenotype					
	MIC (mg/L)	Zone Diameter (mm)			
FOT or TAZ	> 1	< 21 (FOT); < 22 (TAZ)			
MERO	≤ 0.12	≥ 25			
FOX	> 8	< 19			
FOT/CLV and/or TAZ/CLV	NO SYNERGY	NO SYNERGY			
3. ESBL + AmpC-Phenotype					
	MIC (mg/L)	Zone Diameter (mm)			
FOT or TAZ	> 1	< 21 (FOT); < 22 (TAZ)			
MERO	≤ 0.12	≥ 25			
FOX	> 8	< 19			
FOT/CLV and/or TAZ/CLV	SYNERGY	SYNERGY			
			Susceptible		
	MIC (mg/L)	Zone Diameter (mm)			
FOT or TAZ	≤ 1	≥ 21 (FOT); ≥ 22 (TAZ)			
MERO	≤ 0.12	≥ 25			
FOX	≤ 8	≥ 19			

Figure 1: Adapted from EFSA (European Food Safety Authority) and ECDC (European Centre for Disease Prevention and Control), 2020 – The European Union summary report on antimicrobial resistance in zoonotic and indicator bacteria from humans, animals and food in 2017/2018 – and in accordance with the EUCAST guidelines for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance, Version 2.0, July 2017.

The genotype obtained by PCR and/or sequencing may be necessary to correctly categorize a bacterial test strain as either of the categories, ESBL-, AmpC, and/or carbapenemase-producer, but it is not requested as part of this EQA.

Table 2. Interpretive criteria for *Acinetobacter* spp. antimicrobial susceptibility testing

The highlighted antimicrobials could be omitted by the Human Health laboratories.

Antimicrobials	Reference value MIC (µg/mL)			Reference value Disk diffusion (mm)		
	S	I	R	S	I	R
Amikacin, AMK	≤ 16	32	≥ 64	≥ 17	15-16	≤ 14
Cefepime, FEP	≤ 8	16	≥ 32	≥ 18	15-17	≤ 14
Cefotaxime, FOT	≤ 8	16-32	≥ 64	≥ 23	15-22	≤ 14
Ceftazidime, TAZ	≤ 8	16	≥ 32	≥ 18	15-17	≤ 14
Ciprofloxacin, CIP	≤ 1	2	≥ 4	≥ 21	16-20	≤ 15
Colistin, COL	-	≤ 2	≥ 4	NA	NA	NA
Doripenem, DOR	≤ 2	4	≥ 8	≥ 18	15-17	≤ 14
Doxycycline, DOX	≤ 4	8	≥ 16	≥ 13	10-12	≤ 9
Gentamicin, GEN	≤ 4	8	≥ 16	≥ 15	13-14	≤ 12
Imipenem, IMI	≤ 2	4	≥ 8	≥ 22	19-21	≤ 18
Levofloxacin, LEVO	≤ 2	4	≥ 8	≥ 17	14-16	≤ 13
Meropenem, MERO	≤ 2	4	≥ 8	≥ 18	15-17	≤ 14
Minocycline, MIN	≤ 4	8	≥ 16	≥ 16	13-15	≤ 12
Piperacillin/tazobactam, PT4	≤ 16/4	32/4-64/4	≥ 128/4	≥ 21	18-20	≤ 17
Tigecycline, TGC*	≤ 0.5	-	≥ 1	NA	NA	NA
Tobramycin, TOB	≤ 4	8	≥ 16	≥ 15	13-14	≤ 12
Trimethoprim/sulfamethoxazole, SXT	≤ 2/38	-	≥ 4/76	≥ 16	11-15	≤ 10

Reference values are based on *Acinetobacter* spp. breakpoints from CLSI M100, 32nd Ed.

*Reference values are based on *Acinetobacter* spp. clinical breakpoints from www.eucast.org (Tables v. 12.0, 2022).

Table 3. Interpretive criteria for *Staphylococcus aureus* antimicrobial susceptibility testing

The highlighted antimicrobials could be omitted by the Human Health laboratories.

Antimicrobials	Reference value MIC (µg/mL)			Reference value Disk diffusion (mm)		
	S	I	R	S	I	R
Cefoxitin, FOX	≤ 4	-	≥ 8	≥ 22	-	≤ 21
Chloramphenicol, CHL	≤ 8	16	≥ 32	≥ 18	13-17	≤ 12
Ciprofloxacin, CIP	≤ 1	2	≥ 4	≥ 21	16-20	≤ 15
Clindamycin, CLI	≤ 0.5	1-2	≥ 4	≥ 21	15-20	≤ 14
Erythromycin, ERY	≤ 0.5	1-4	≥ 8	≥ 23	14-22	≤ 13
Fusidate, FUS*	≤ 1	-	≥ 2	≥ 24	-	≤ 23
Gentamicin, GEN	≤ 4	8	≥ 16	≥ 15	13-14	≤ 12
Kanamycin, KAN*	≤ 8	-	≥ 16	≥ 18	-	≤ 17
Linezolid, LZD	≤ 4	-	≥ 8	≥ 21	-	≤ 20
Penicillin, PEN	≤ 0.12	-	≥ 0.25	≥ 29	-	≤ 28
Quinupristin/dalfopristin, SYN	≤ 1	2	≥ 4	≥ 19	16-18	≤ 15
Rifampin, RIF	≤ 1	2	≥ 4	≥ 20	17-19	≤ 16
Sulfamethoxazole, SMX	≤ 256	-	≥ 512	≥ 17	13-16	≤ 12
Tetracycline, TET	≤ 4	8	≥ 16	≥ 19	15-18	≤ 14
Trimethoprim, TMP	≤ 8	-	≥ 16	≥ 16	11-15	≤ 10
Vancomycin, VAN	≤ 2	4-8	≥ 16	NA	NA	NA

Reference values are based on *Staphylococcus aureus* breakpoints from CLSI M100, 32nd Ed.

*Reference values are based on *Staphylococcus aureus* clinical breakpoints from www.eucast.org (Tables v. 12.0, 2022).

4 REPORTING OF RESULTS AND EVALUATION

We recommend that you write your results in the enclosed test forms and that you read carefully the description in paragraph 5 before entering your results in the Informatics Module. If the same reference strain is used for different pathogens, please enter the results (even if the same) for all the pathogens. The Informatics Module will allow you to view and print a report with your reported results. The scores for the results will be released after the result submission deadline where you will be able to access the evaluation of your results. Results in agreement with the expected interpretation are categorised as ‘4’ (correct), while results deviating from the expected interpretation are categorised as ‘3’ (incorrect, minor), ‘1’ (incorrect, major) or ‘0’ (incorrect, very major).

Results must be submitted no later than June 10th 2022.

If you have trouble in entering your results, please contact the EQA Coordinator directly, explaining the issues that you encountered:

Patrícia T. dos Santos
National Food Institute, Technical University of Denmark
Kemitorvet, Building 204, DK-2800 Lyngby – DENMARK
E-mail: pado@food.dtu.dk

Direct communication with the EQA Coordinator must be in English.

5 HOW TO SUBMIT RESULTS VIA THE INFORMATICS MODULE

The ‘Guideline for reporting results in the EQAsia Informatics Module’ is available for download directly from the [EQAsia website](#). Please follow the guideline carefully.

Access the Informatics Module (incognito window) using [this address](#). See below how to login to the Informatics Module.

When you submit your results, remember to have by your side the completed test forms (template available for download from the [EQAsia website](#)).

Do not hesitate to contact us if you have trouble with the Informatics Module.

Before finally submitting your input for all the organisms, please ensure that you have filled in all the relevant fields as **you can only ‘finally submit’ once!** ‘Final submit’ blocks data entry.

Login to the Informatics Module:

When first given access to login to the Informatics Module, your **personal loginID and password** is sent to you by email.

Note that the primary contact person for a participating institution is registered both as primary and secondary contact. Should you like to add another person as the secondary contact, please contact pado@food.dtu.dk

--- --- ---